

Continental SCIENTIFIC

Extra Wide Single Section General Purpose Laboratory/Pharmacy Refrigerators

Utilizing Environmentally Safe, CFC-free, R-134a Refrigerant



S1RX-6D



S1RX

Shown w/ optional
Stainless Steel Drawers

Standard Features

- Modern, State-of-the Art Styling
- Performance Rated Refrigeration System Utilizing Environmentally Safe R-134a Refrigerant
- 3" Non-CFC Polyurethane Foam Insulation
- Chrome-Plated Flow Line Handles
- Cam-Action Lift-Off Hinges
- Self-Closing Doors
- Automatic Interior Lighting
- Magnetic Snap-In Door Gaskets
- Cylinder Lock In Each Door
- Automatic Electric Condensate Evaporator
- 5" Casters
- Heavy-Duty, Epoxy-Coated Steel Shelves
- External Digital Thermometer
- Adjustable Temperature Control, Range +2°C to +8°C, Factory Present at +4°C


Standard (S/S Front, Aluminum End Panels & Interior)
Suffix SA (S/S Exterior, Aluminum Interior)
Suffix SS (S/S Exterior & Interior)

Optional Features & Accessories

- Legs
- Pass-Thru Doors
- Additional Shelves
- Rehinging of Doors
- Half Height Doors
- Stainless Steel Drawers
- 2" Access Port(s)
- 220/50/1 Voltage
- 220/60/1 Voltage
- Chromatography Pole Assembly
- Remote
- 7 Day Temperature Chart Recorder
- Audible/Visual Temperature Alarm
- Export Crating

Project Name: _____




Item Number: _____



Extra Wide Single Section General Purpose Laboratory/Pharmacy Refrigerators

SPECIFICATIONS

CHARACTERISTICS DIMENSIONAL DATA	<u>Model</u> S1RX One Section		<u>Model</u> S1RX-GD One Section	
	Reach-In	Pass-Thru	Reach-In	Pass-Thru
Net Capacity (cu. ft.)	30	30	30	30
Width, Overall (in.)	36 ¹ / ₄	36 ¹ / ₄	36 ¹ / ₄	36 ¹ / ₄
Depth, Overall (Incl. Hdls.) (in.)	35 ¹ / ₂	39	35 ¹ / ₂	39
Depth (Less Door(s)) (in.)	32	32	32	32
Depth (Door(s) Open 90°) (in.)	66	100	66	100
Clear Door Width (in.)	29 ⁵ / ₈	29 ⁵ / ₈	29 ⁵ / ₈	29 ⁵ / ₈
Clear Full-Door Height (in.)	58 ⁵ / ₈	58 ⁵ / ₈	58 ⁵ / ₈	58 ⁵ / ₈
Clear Half-Door Height (in.)	27 ¹ / ₂	27 ¹ / ₂	27 ¹ / ₂	27 ¹ / ₂
Height, Overall (inc. 5" Casters) (in.)	82 ¹ / ₄	82 ¹ / ₄	82 ¹ / ₄	82 ¹ / ₄
No. of Full Door(s)	1	2	1	2
No. of Half Door(s)	2	4	2	4
No. of Shelves	3	3	4	4
Shelf Area (sq. ft.)	24.4	24.4	24.4	24.4
REFRIGERATION DATA				
Condensing Unit Size (H.P.)	¹ / ₄	¹ / ₄	¹ / ₃	¹ / ₃
Capacity (BTU/Hr)*	2190	2190	2810	2810
ELECTRICAL DATA				
Cabinet Specifications	115/60/1	115/60/1	115/60/1	115/60/1
Feed Wires (Incl. Ground)	3	3	3	3
Condensing Unit Specs.	115/60/1	115/60/1	115/60/1	115/60/1
Total Amperes	7.1	7.4	11.0	11.3
Cabinet Amperes	1.8	2.1	3.8	4.1
Defrost Amperes	-	-	-	-
Power Supply Cord	yes	yes	yes	yes
SHIPPING DATA				
Height - crated (in.)	82	82	82	82
Width - crated (in.)	42	42	42	42
Depth - crated (in.)	42	42	42	42
Volume - crated (cu.ft.)	87	87	87	87
Weight Std - crated (lbs.)	430	435	430	435
Weight SS - crated (lbs.)	490	495	490	495
Weight Std - uncrated (lbs.)	330	335	330	335
Weight SS - uncrated (lbs.)	390	395	390	395



* Rating @ +25°F Evaporator, 90° Ambient

Note: All specifications subject to change without notice.

Note: If the cabinet is located directly against wall and/or under low ceiling, a minimum clearance of 2" is required.

Continental

SCIENTIFIC

A Division of National Refrigeration & A/C Products
539 Dunksferry Road • Bensalem, PA 19020
(215) 244-1400 • (800) 523-7138 • Fax: (215) 244-9579
www.continentalscientific.com