

Continental *SCIENTIFIC*

Sliding Glass Door General Purpose Laboratory/Pharmacy Refrigerators

Utilizing Environmentally Safe, CFC-free, R-134a Refrigerant



Standard Features

- Modern, State-of-the Art Styling
- Performance Rated Refrigeration System Utilizing Environmentally Safe R-134a Refrigerant
- 3" Non-CFC Polyurethane Foam Insulation
- Chrome-Plated Flow Line Handles
- Double-Pane, Tempered Sliding Glass Doors w/ Lock
- Self-Closing Doors
- Fluorescent Interior Lighting
- Magnetic Snap-In Door Gaskets
- Automatic Electric Condensate Evaporator
- 5" Casters
- Heavy-Duty, Epoxy-Coated Steel Shelves
- External Digital Thermometer
- Adjustable Temperature Control, Range +2°C to +8°C, Factory Present at +4°C


Standard (S/S Front, Aluminum End Panels & Interior)
Suffix SA (S/S Exterior, Aluminum Interior)
Suffix SS (S/S Exterior & Interior)

Optional Features & Accessories

- Legs
- Pass-Thru Doors
- Additional Shelves
- Rehinging of Doors
- Half Height Doors
- Audible/Visual Temperature Alarm
- 220/50/1 Voltage
- 220/60/1 Voltage
- Chromatography Pole Assembly
- Remote
- 7 Day Temperature Chart Recorder
- 2" Access Port(s)
- Vapor Proof Duplex Outlet
- Pass-Thru Doors
- Export Crating

Project Name: _____



Item Number: _____



Sliding Glass Door General Purpose Laboratory/Pharmacy Refrigerators

SPECIFICATIONS

	<u>Model</u> S2R-SGD	<u>Model</u> S2R-SGD-PT
	Reach-In	Pass-Thru
CHARACTERISTICS		
DIMENSIONAL DATA		
Net Capacity (cu. ft.)	48	48
Width, Overall (in.)	52	52
Depth, Overall (Incl. Hdl.s) (in.)	34	36
Depth (Less Door(s)) (in.)	32	32
Clear Door Width (in.)	19	19
Clear Full-Door Height (in.)	58 ¹ / ₄	58 ¹ / ₄
Clear Half-Door Height (in.)	27	27
Height, Overall (inc. 5" Casters) (in.)	82 ¹ / ₄	82 ¹ / ₄
No. of Full Door(s)	2	4
No. of Half Door(s)	4	8
No. of Shelves	8	8
Shelf Area (sq. ft.)	40.8	40.8
REFRIGERATION DATA		
Condensing Unit Size (H.P.)	1/2	1/2
Capacity (BTU/Hr)*	4010	4010
ELECTRICAL DATA		
Cabinet Specifications	115/60/1	115/60/1
Feed Wires (Incl. Ground)	3	3
Condensing Unit Specs.	115/60/1	115/60/1
Total Amperes	13.6	14.3
Cabinet Amperes	4.1	4.8
Defrost Amperes	-	-
Power Supply Cord	yes	yes
SHIPPING DATA		
Height - crated (in.)	82	82
Width - crated (in.)	55	55
Depth - crated (in.)	37	37
Volume - crated (cu.ft.)	97	97
Weight Std - crated (lbs.)	530	610
Weight SS - crated (lbs.)	610	680
Weight Std - uncrated (lbs.)	360	430
Weight SS - uncrated (lbs.)	440	510



* Rating @ +25°F Evaporator, 90° Ambient

Note: All specifications subject to change without notice.

Note: If the cabinet is located directly against wall and/or under low ceiling, a minimum clearance of 2" is required.

Continental
SCIENTIFIC

A Division of National Refrigeration & A/C Products
539 Dunksferry Road • Bensalem, PA 19020
(215) 244-1400 • (800) 523-7138 • Fax: (215) 244-9579
www.continentalscientific.com